



Be Right™



## EDTA Digital Titrator Cartridge, 0.714 M

**Termékszám:**

**1495901**

**Veszélyes**

HUF Ár (nettó):

Elérhetőségünk

Elérhető



Ready-to-use titrant solution for the analysis of Total Hardness, 10-100+ G.d.H. CaCO<sub>3</sub>.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Műszaki adatok

Méréstartomány:	10 - 200 G,d,h,
Módszer név:	EDTA (CalVer2)
Módszer száma:	8204
Paraméter:	Hardness, Total
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titrló:	0.714 M EDTA