



Iodate-Iodide Digital Titrator Cartridge, 0.3998 N

Termékszám:

1496101

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



Ready-to-use titrant solution for the analysis of Sulfite, 4 - >400mg/L SO₃.

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	0 - 800 mg/L SO ₃ ²⁻
Módszer név:	Iodate-Iodide
Paraméter:	Sulfite
Platform:	Digitális titrátor
Titrló:	0.3998 N KIO ₃ -KI