



## Maintenance kit for 5500sc Silica Analyser, 6 channels

Termékszám:	6788303
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### Műszaki adatok

A doboz tartalma:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions
Product Kit:	5500 sc Silica Analyser, 6 Channels

---

### A doboz tartalma

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions