



Iron color disc test kit, model IR-18B

Termékszám:

146401

Veszélyes

HUF Ár (nettó):
Elérhető

Elérhetőségünk



Color Disc Kits - Single Parameter

Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	NA
Méréstartomány (2):	NA
Módszer név:	Color disc/ 1,10 Phenanthroline
Paraméter:	Iron, medium range - as Fe
Platform:	Color Disc
Range:	0.25 - 7 mg/L Fe
Smallest increment:	0.25
Szállítási tömeg:	1
Szünet mód:	NA
Tartó stílus:	B
Teszttek száma:	100
Tesztkészlet model:	IR-18B
Típus:	IR-18B