



Be Right™



## Test Kit, chlorine (free & total) & pH, model CN-67

**Termékszám:**

**1411100**

**Veszélyes**

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

**One of Hach's most popular test kits**

**Includes color discs and a comparator for monitoring chlorine and pH**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Műszaki adatok

Alkalmazás:	Pool and Spa
Méréstartomány:	6.6 - 8.4 pH units
Méréstartomány (2):	0 - 3,4 mg/L
Platform:	Color Disc
Tesztek száma:	300