



Hardness, iron, and pH test kit, model HA-62A

Termékszám:

183701

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Veszélyes



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Alkalmazás:	Water Quality
Méréstartomány:	6.6 - 8.4 pH units
Paraméter:	Hardness, Iron, and pH (6.6-8.4)
Platform:	Color Disc
Tesztek száma:	100
Típus:	HA-62A