



Test kit, manganese color disc, model MN-PAN

Termékszám:

2350800

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	0 - 0.70 mg/L Mn
Módszer név:	Color disc/PAN
Paraméter:	Manganese - As Mn
Parameter:	Manganese
Platform:	Color Disc
Smallest increment:	0.05
Standard mérési tartomány:	0 - 0.70 mg/L Mn
Tartó stílus:	D
Tesztek száma:	50
Tesztkészlet model:	MN-PAN
Típus:	MN-PAN