



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Termékszám:

143601

Veszélyes

HUF Ár (nettó):
Elérhető

Elérhetőségünk



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Méréstartomány:	2 - 40 mg/L
Méréstartomány (2):	1.25 - 25 mg/L
Módszer név:	Drop count titration/ Sodium hydroxide Phenolphthalein
Paraméter:	Carbon Dioxide - As CO ₂
Platform:	Drop Count Titration
Smallest increment:	1.25 2 5
Tartó stílus:	D
Tesztek száma:	200
Tesztkészlet model:	CA-23