



Be Right™



Test kit cyanide, model CYN-3

Termékszám:

201002

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Veszélyes



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Műszaki adatok

Méréstartomány:	0 - 0.3 mg/L CN ⁻
Módszer név:	Color disc/ Pyridine-Pyrazolone
Paraméter:	Cyanide - As free CN ⁻
Parameter:	Cyanide
Platform:	Color Disc
Smallest increment:	0.02
Standard mérési tartomány:	0 - 0.3 mg/L CN ⁻
Tartó stílus:	R
Teszttek száma:	100
Tesztkészlet model:	CYN-3
Típus:	CYN-3