



Total Hardness Test Kit, Model 5-EP

Termékszám:

145400

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model 5-EP test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizes drop count titration with a range of 1-30 gpg. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	1 - 30 CaCO ₃
Módszer név:	Drop count titration/ EDTA-powder pillows
Paraméter:	Hardness, Total - As CaCO ₃
Platform:	Csepp számlálós titrátor
Smallest increment:	1 gpg
Tartó stílus:	B
Tesztek száma:	100
Tesztkészlet model:	5-EP