



Be Right™



Test kit, alkalinity, model AL-DT

Termékszám: 2063700
HUF Ár (nettó): Elérhetőségünk
Szállítás 2 héten belül

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Legkisebb beosztás:	0.1-10
Méréstartomány:	10 - 4,000 mg/L CaCO ₃
Módszer név:	Digital Titrator/Sulfuric acid
Paraméter:	Alkalinity - (P) & Total (MO) as CaCO ₃
Platform:	Digitális titrátor
Szállítási tömeg:	3.5
Tartó stílus:	F
Tesztek száma:	100
Tesztkészlet model:	AL-DT