



Be Right™



Alkalinity test kit AL-AP, 5-100/20-400 mg/L CaCO₃ 100 pcs

Termékszám: 2444301
HUF Ár (nettó): Elérhetőségünk
Szállítás 1 héten belül

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Méréstartomány:	5 - 100 mg/L
Méréstartomány (2):	20 - 400 mg/L
Módszer név:	Drop count titration/Sulfuric acid
Paraméter:	Alkalinity - (P) & Total (MO) as CaCO ₃
Platform:	Csepp számlálós titrátor
Smallest increment:	5 20
Szállítási tömeg:	1
Tartó stílus:	B
Tesztek száma:	100
Tesztkészlet model:	AL-AP, MG-L