



POCKET Colorimeter II Colorimeter Test Kit for Manganese (HR) analysis

Termékszám: 5870015

Veszélyes

ELAVULT CIKK
Ez a termék már nem elérhető.



Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Manganese, 0.20-20.0 mg/L Mn, Periodate Oxidation Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Műszaki adatok

Borítás besorolás:	IP67 (elemtartó nélkül)
Detektor:	Silicon detector
Hullámhossz:	Fixed wavelength ± 2 nm varies with model
Kijelző:	LCD, backlit
Megfelelőségi igazolások:	CE
Méréstartomány:	0.2 - 20.0 mg/L Mn
Méretek (MxSZxM):	155 mm x 61 mm x 35 mm
Módszer név:	Cold periodate
Módszer száma:	8034
Paraméter:	Manganese
Smallest increment:	0.1
Súly:	0,23 kg
Tartó stílus:	D
Tesztek száma:	100