



Be Right™



Sodium Hydroxide Digital Titrator Cartridge, 0.3636 N

Termékszám:

1437801

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	10 - 100 mg/L CO ₂
Módszer név:	Phenolphthalein
Módszer száma:	8205
Paraméter:	Carbon Dioxide
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titráló:	0.3636 N NaOH