



Alkalinity Reagent Set for Digital Titrator

Termékszám:

2271900

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Ready to use titrants in easily exchangeable, resealable cartridges. For reliable results with the Digital Titrator in the field and in the laboratory.

For alkalinity determination by sulfuric acid titration with the Digital Titrator. Hach method 8203. Range: 10-4,000 mg/L CaCO₃ Set contains Phenolphthalein Indicator Powder Pillows, Bromcresol Green-Methyl Red Indicator Powder Pillows and Sulfuric Acid Digital Titrator Cartridges, 0.1600 N and 1.600 N. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Különleges jegyzet:

Set contains: Bromcresol Green Methyl Red Indicator Powder Pillows (94399), Phenolphthalein Indicator (94299), Sulfuric Acid Titration Cartridge 1.600 N (1438901), Sulfuric Acid Titration Cartridge 0.1600 N (1438801)

Leírás:

Alkalinity Reagent Set for Digital Titrator

Méréstartomány:

10 - 4000 mg/L CaCO₃

Módszer név:

Phenolphthalein and Total Alkalinity

Módszer száma:

8203

Paraméter:

Alkalinity

Platform:

Digitális titrátor

Tesztek száma:

100