



Stablcal Stabilizált formazin zavarosság standard 1.0 NTU

Termékszám:

2659842

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Használatra kész, stabilizált formazin elsődleges standard oldat, 100 mL, 1,0 NTU.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Műszaki adatok

Koncentráció:	1 NTU
Platform:	Oldat
Térfogat/ csomagméret:	ml