



## Stablcal Turbidity Standard, 800 NTU, 100 mL

**Termékszám:**

**2660542**

**Veszélyes**

HUF Ár (nettó):  
Elérhető

Elérhetőségünk



Ready to use, stabilised Formazine Primary standard solution in 100 mL bottle, 800 NTU.

USEPA accepted for calibration purposes. Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. For TL23, TU52, 2100Q, 2100N, 2100AN series.

### Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Műszaki adatok

Koncentráció:	800 NTU
Platform:	Oldat
Térfogat/ csomagméret:	100 mL