



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 500 mL

Termékszám:

2947549

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 500 mL bottle

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Műszaki adatok

EPA compliant:	Igen
Méréstartomány:	0.02 - 2.00 mg/L F ⁻
Módszer név:	SPADNS 2
Módszer száma:	10225
Platform:	Oldat
Tesztek száma:	125