



Be Right™



EDTA Digital Titrator Cartridge, 0.1428 M

Termékszám:

1496001

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Ready-to-use titrant solution for the analysis of Total Hardness, 1-16 G.d.H. CaCO₃.

For the determination of Calcium and Total Hardness. 0.1428 M. Digits x 0.1785 = mg/L CaCO₃ in 100 mL, Digits x 0.01 = Gdh.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	1 - 16 G,d,h,
Módszer név:	EDTA (CalVer® 2)
Paraméter:	Hardness, Total
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titráló:	0.1428 M EDTA