



Be Right™



## Filter membrane, EZ-Pak, 0.22 micron, 47 mm, white, gridded, 600/pk (millipore)

**Termékszám:** 2916900  
**HUF Ár (nettó):** Elérhetőségünk  
**Szállítás 2 héten belül**

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

---

### Műszaki adatok

|                        |                       |
|------------------------|-----------------------|
| Anyag:                 | Mixed-Cellulose Ester |
| Átmérő:                | 47 mm mm              |
| Szemcseméret:          | 0.22 µm               |
| Térfogat/ csomagméret: | 600                   |