



## Filter kit, field, coliform & E.coli

**Termékszám:**

**2586400**

**Veszélyes**

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



### MEL/m-ColiBlue24® Field Filtration Lab

The Coliform and E. coli Field Filtration Kit features Hach's m-ColiBlue24® Broth and Millipore's Portable filter assembly, Push Fit funnels and Millipore membranes to simplify membrane filtration testing in the field. m-ColiBlue24® Broth simultaneously detects total coliforms and E. coli in 24 hours.

The kit includes all the necessary materials and apparatus for convenient, efficient field operations, and comes in a carrying case for easy storage.

Contains: 40 m-ColiBlue24® broth ampules, 48 Push Fit funnels and S-Pak filters (Millipore), 1 Microfil Stand Alone support (Millipore), 140 cc syringe with tubing, 50 petri dishes (50 mm, sterile) with nutrient pad, 50 sterile, individually wrapped 10 ml pipets (poly), 50 sterile sample bags with dechlorinating agent, alcohol burner and 100 ml of isopropyl alcohol, illuminated magnifier (2.5 and 5x), instructions for procedure.

### Consumables and media included

### EPA-approved field method

---

### Műszaki adatok

|                |        |
|----------------|--------|
| Érzékenység:   | 1      |
| Tesztek száma: | Varies |