



POCKET Colorimeter II Koloriméter készlet formaldehid küvettatesztekhez

Termékszám: 5953000V.08

ELAVULT CIKK
Ez a termék már nem elérhető.

Portable Colorimeter programmed for the determination of Formaldehyde with cuvette tests. In carrying case, with batteries and operating instructions.

Following Hach cuvette test can be evaluated with this instrument:

LCK325, 0.5-10.0 mg/L H₂CO

Cuvette test has to be ordered separately.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to lightning and large figures on the display

Reliable results with no mains connection

Rugged construction yet lightweight

Műszaki adatok

Méréstartomány: 0.5 - 10.0 mg/L H₂CO

Paraméter: Formaldehid (az LCK325 küvettához alkalmas)