



Beaker, polypropylene, low form, 50 mL

Termékszám:	108041
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Műszaki adatok

Anyag:	Polypropylene, Low Form
Anyaga:	Polypropylene
Átmérő, külső:	51.1 mm with spout
Egység:	1 each
Kapacitás:	50 mL
Térfogat/ csomagméret:	1 each