



Be Right™



Funnel, analytical, 254 mL approx. Volume

Termékszám:	108370
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

Műszaki adatok

Térfogat/ csomagméret:	1 each
------------------------	--------