



Flask, volumetric, polypropylene, 250 mL

Termékszám: 1406046

ELAVULT CIKK

Ez a termék már nem elérhető.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Műszaki adatok

Anyaga:	Polypropylene
Egység:	1 each
Kapacitás:	250 mL
Tanúsítványok:	Not Applicable