



Standard solution, phosphate, 500 mg/L as PO₄ (NIST), pk/16 - 10-mL voluette[®] ampules

Termékszám:	1424210
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

For accuracy checks in phosphate (phosphorus) determinations

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Műszaki adatok

Koncentráció:	500 mg/L
Leírás:	Phosphate Standard Solution
Platform:	Oldat
Térfogat/ csomagméret:	mL/16