



Be Right™



m-FC broth, dehydrated, 100 g

Termékszám:

1462426

ELAVULT CIKK

Ez a termék már nem elérhető.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Műszaki adatok

Media Type: m-FC Broth

Térfogat/ csomagméret: 100 g