



Cyanuric acid test kit, model CY-3

Termékszám: 185102
HUF Ár (nettó): Elérhetőségünk
Szállítási dátum nem látható

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Alkalmazás:	Pool and Spa
Méréstartomány:	20 - 100 mg/L
Módszer név:	Turbidimetry
Paraméter:	Cyanuric Acid
Smallest increment:	20
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	CY-3
Típus:	CY-3