



Be Right™



## Digital titrator cartridge, phenylarsine oxide (PAO), 0.00564 N

Termékszám:

199901

Veszélyes

### ELAVULT CIKK

Ez a termék már nem elérhető.



For low level chlorine determination by amperometric titration, using the Digital Titrator

#### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

#### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

#### Versatile

More than 50 parameter/range combinations are available

#### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

## Műszaki adatok

EPA compliant:	Igen
Méréstartomány:	0 - 1000 µg/L Cl <sub>2</sub>
Méréstartomány (2):	0 - 1000 µg/L Cl <sub>2</sub>
Módszer név:	Amperometric, Forward
Módszer száma:	10026
Paraméter:	Chlorine, Free
Platform:	Digitális titrátor
Titráló:	0.00564 N PAO