



Test kit, hardness & iron color disc, model HA-77

Termékszám:

202300

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Veszélyes



Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

The kit to use for more resolution in a visual test

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Alkalmazás:	Water Quality
Paraméter:	Hardness and Iron
Platform:	Color Disc
Típus:	HA-77