



Test kit, general boiler water, drop count titration, model AL-94 MG-L

Termékszám:

222403

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Veszélyes



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Műszaki adatok

Alkalmazás:	Boiler Water
Paraméter:	Drop Count Titration
Platform:	Color Disc
Típus:	AL-94 MG-L