



Test kit, water conditioning, model SO-1

Termékszám: 225300

ELAVULT CIKK

Ez a termék már nem elérhető.

Hardness/Soap "Demo" for Sales Aid in Water Conditioning. Use drop-count methods to demonstrate the relationship of hardness to soap consumption. Note: While it provides a quick hardness/soap check, it is less accurate than the hardness titration supplied by Model 15-E (PN 148201) and Model AD-16E (PN 148410) water conditioning test kits.

Uses drop-count methods to demonstrate the relationship of hardness to soap consumption

Provides an easy, economical, and fast hardness/soap check

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	1 - 20 gpg
Paraméter:	Soap Demonstration
Platform:	Drop Count Titration