



Test kit, chlorine (total), model CN-65

Termékszám:

225401

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Legkisebb beosztás:	1
Méréstartomány:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Módszer név:	Drop count titration/ Thiosulfate
Paraméter:	Chlorine, mid and high range - As Total Cl ₂
Platform:	Drop Count Titration
Smallest increment:	0.2
	1
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	CN-65