



Digital titrator cartridge, sodium thiosulfate, 0.2068N

Termékszám:	2267601
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

For determination of Chromate by Digital Titrator titration. 0.2068 N Na₂S₂O₃. Digits x 0.2 = mg/L CrO₄ in a 50 mL sample.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	20 - 400 mg/L CrO ₄
Módszer név:	Iodometric
Módszer száma:	8211
Paraméter:	Chromate
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titrló:	0.2068 N Na ₂ S ₂ O ₃