



Reagent set, carbon dioxide, digital titrator

Termékszám:

2272700

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



For determination of Carbon Dioxide by Digital Titrator titration. Hach Method 8205. Range: 10 - 1000 mg/L as CO₂ Set includes: Phenolphthalein Pillows and two Sodium Hydroxide Digital Titrator Cartridges (0.3636 N and 3.3636 N). Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	10 - 1000 mg/L CO ₂
Módszer száma:	8205
Platform:	Digitális titrátor
Tesztek száma:	40 - 100