



Be Right™



Test kit, ascorbic acid, model ASC-1

Termékszám:

2308100

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Méréstartomány:	10 - 200 mg/L
Módszer név:	Drop count titration/ Iodometric
Paraméter:	Ascorbic Acid
Platform:	Drop Count Titration
Smallest increment:	10
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	ASC-1