



## Butterfield's buffered phosphate, 90 mL PK/25

<b>Termékszám:</b>	<b>2319125</b>
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Butterfields Buffered Phosphate. 90 mL bottle. pk/25.

### **Complete Solutions for Testing Recreational Water**

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

### **Complete Solutions for Testing Wastewater**

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

### **Complete Solutions for Testing Drinking Water**

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

---

## **Műszaki adatok**

Térfogat:	90 mL
Térfogat/ csomagméret:	25/pk