



Test kit, molybdate, model MO-LR

Termékszám:

2359300

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Műszaki adatok

Méréstartomány:	Single Range
Módszer név:	Color disc/
	Ternary complex
Paraméter:	Molybdate - As Molybdenum
Parameter:	Molybdate
Platform:	Color Disc
Smallest increment:	0.25
Standard mérési tartomány:	0 - 3.0 mg/L Mo
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	MO-LR
Típus:	MO-LR