



Agar LES, m-Endo, dehydrated 500 g

Termékszám:

2406634

Veszélyes

ELAVULT CIKK

Ez a termék már nem elérhető.



Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Műszaki adatok

Media Type:	m-Endo LES Agar
Térfogat/ csomagméret:	500 g