



Test kit, nitrogen, ammonia, model NI-SA

Termékszám:

2428700

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Műszaki adatok

Méréstartomány:	0 - 2.4 mg/L NH ₃ -N
Módszer név:	Color disc/Salicylate
Paraméter:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Platform:	Color Disc
Smallest increment:	0.1
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	NI - SA