



pH Paper, 0 - 14 pH Range, 100 tests

Termékszám:	2601300
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples. Measures 0 - 14 pH.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Műszaki adatok

Egység:	100 tests
Legkisebb beosztás:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Méréstartomány:	0 - 14 pH units
Módszer száma:	10056
Paraméter:	pH
Platform:	Tesz csík
Tesztek száma:	100