



Cylinder, graduated, 1000 mL, w/ handles

Termékszám:	2612953
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered ”double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Műszaki adatok

Anyaga:	Polypropylene
Kapacitás:	1000 mL
Mérés módszer:	To Deliver/To Contain
Tanúsítványok:	None
Térfogat/ csomagméret:	1 each
Típus:	Dupla metrikus