



## Bottle, wash, diamond alert, 500 mL, Variety 6/pk

|                         |                |
|-------------------------|----------------|
| Termékszám:             | 2613049        |
| HUF Ár (nettó):         | Elérhetőségünk |
| Szállítás 2 héten belül |                |

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

---

### Műszaki adatok

|                        |               |
|------------------------|---------------|
| Anyag:                 | LDPE          |
| Egység:                | 6 pk          |
| Kapacitás:             | 500 mL        |
| Stem:                  | Blue Cap/Stem |
| Szín:                  | Transparent   |
| Térfogat/ csomagméret: | 6 pk          |