



Flask, volumetric, class A, certified, 1000 mL

Termékszám: 2636653

ELAVULT CIKK

Ez a termék már nem elérhető.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

Műszaki adatok

Anyaga:	Class A Glass
Egység:	1 each
Kapacitás:	1000 mL
Tanúsítványok:	Yes