



Clamp, Extension Thermometer

Termékszám:	2656200
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Replacement clamp for use with the Cold Vapor Mercury Test. Hold thermometers and other glass apparatus firmly with this lightweight, safety-adjust clamp. Holds 6.4 - 12.7 mm (1/4 - 1/2") diameter glass tubing/thermometers up to 178 mm (7") from support rod. Tightening wing nut applies tension to the nickel-plated phosphor bronze jaws. Overall length: 210 mm (8 1/4"). If attaching to rods, a 90° clamp holder is required. 1/pk