



Reagent set monochloramine, (indophenol method) HR TNT

Termékszám:

2805145

Veszélyes

ELAVULT CIKK

Ez a termék már nem elérhető.



Javasolt cserék:

- 2802246

For determination of Monochloramine by the Indophenol Test 'N Tube[™] method. Hach Method 10172. Range: 0.1 - 10.0 mg/L as Cl₂. Sample Size: 2 mL. Up to 50 tests/reagent set.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Műszaki adatok

Lábjegyzet:	+ Dependent on the number of performed reagent blanks.
Méréstartomány:	0.1 - 10.0 mg/L
Módszer név:	Indophenol
Módszer száma:	10172
Műszer:	All except DR 900, PC II
Paraméter:	Chloramine, Mono
Platform:	Test 'N Tube
Tesztek száma:	50