



Vacuum filter microfunnel, pk/50.

Termékszám:	2831500
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metricel with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Műszaki adatok

Anyag:	Mixed-Cellulose Ester
Átmérő:	47 mm
Szemcseméret:	0.45 µm
Térfogat/ csomagméret:	50