



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Termékszám:	2939103
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Műszaki adatok

Anyag:	PP
Body Material:	Plastic
Kapacitás:	125 mL
Package size UOM:	12pk
Térfogat/ csomagméret:	12 /pk