



Fixed Volume, Micropipette, 5 µL

Termékszám:	2946001
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and colour-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Műszaki adatok

Kapacitás:	5 µL
Pontosság:	1,5 %
Térfogat/ csomagméret:	1 each