



Fixed volume, Micropipette, 1000 µL

Termékszám:	2946010
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Műszaki adatok

Kapacitás:	1000 µL
Pontosság:	0,5 %
Térfogat/ csomagméret:	1 each